

CP	BULL	MASCU	THICK	MUSCLE	RATING	DISP	COMMENTS
PAD	1106X	3.1	3.1	3.1	2.9	4.9	big spot left eye
PMS	1145X	4.1	4.1	4.0	4.0	5.0	
PMS	1146X	4.1	4.0	4.0	4.0	5.0	
PMS	1147X	4.9	4.1	4.0	4.1	5.0	
PMS	1149X	4.0	3.9	3.9	3.2	4.0	
PMS	1151X	4.0	4.0	4.0	4.0	4.9	
PMS	1152X	4.9	4.9	4.9	4.9	5.0	? Scrotum
PMS	1153X	3.1	3.0				left rear steps out; CULL
PMS	1154X	3.1	4.1	4.0	4.1	5.0	
PMS	1155X	3.0	3.1	3.0	2.9	4.0	
PMS	1156X	3.0	3.9	3.9	3.9	4.9	
PMS	1157X	3.1	3.9	3.9	3.0	4.1	ears dead dogs
OLC	2207X	4.1	4.1	4.1	4.1	2.9	
OLC	2209X	4.0	4.0	4.0	4.0	4.1	
D-J	2701X	5.0	5.0	4.9	5.1	5.0	
D-J	2702X	4.1	4.9	4.9	5.0	4.0	
D-J	2704X	4.9	4.2	5.0	4.2	4.1	
D-J	2705X	4.0	4.0	4.9	4.2	4.1	spot left eye; tall
D-J	2751X	4.0	4.2	4.1	4.2	3.0	
D-J	2752X	3.1	3.9	4.0	3.9	4.9	tall
D-J	2754X	3.1	4.1	5.0	4.9	4.0	
D-J	2755X	4.1	4.1	4.0	4.1	3.9	
D-J	2756X	4.0	5.0	5.0	5.0	4.1	
D-J	2757X	4.9	4.9	4.9	4.1	4.0	
D-J	2758X	3.1	3.9	4.0	3.2	5.0	
D-J	2759X	5.0	5.0	5.0	5.1	5.0	
D-J	2760X	2.9	2.9	3.1	2.8	4.0	

CP	BULL	MASCU	THICK	MUSCLE	RATING	DISP	COMMENTS
D-J	2763X	4.1	4.1	4.9	4.1	5.0	
LLL	5903X	4.1	4.0	4.0	4.1	5.0	
LLL	5906X	5.0	5.0	5.0	5.0	4.1	
LLL	7300X	3.0	2.9	3.9	2.8	4.1	brockle face
R2R	7754X	3.1	3.0	3.0	3.1	4.0	#2
R2R	7755X	4.0	3.9	2.9	4.0	3.9	#1
R2R	7761X	4.0	2.1	2.1	2.1	3.0	#3
34							

CP	BULL	MASCU	THICK	MUSCLE	RATING	DISP	COMMENTS
PAD	4065X	3.1	3.1	4.0	3.9	4.1	
PAD	4066X	3.0	3.9	3.9	2.9	4.0	
PWW	4110X	4.0	4.9	4.0	4.0	4.9	short
PWW	4111X	4.1	4.1	4.2	4.0	4.0	high tail head
PWW	4112X	5.0	4.9	5.0	5.1	4.0	
PWW	4113X	5.1	4.9	4.9	5.0	4.1	
PWW	4114X	5.0	4.9	4.9	4.9	3.0	little man syndrome
PWW	4115X	4.9	4.1	4.1	4.2	4.1	
PWW	4116X	5.0	4.9	4.0	4.9	4.0	
PWW	4117X	5.1	4.1	4.1	4.1	4.9	tail head high
PWW	4118X	5.0	4.1	4.0	4.0	4.0	
PWW	4119X	5.0	4.0	4.0	4.0	3.0	
PWW	4120X	5.1	4.9	4.1	5.0	4.0	
WFF	4503X	5.0	4.9	4.9	5.0	5.0	
WFF	4504X	4.9	4.8	4.0	4.8	5.0	
WFF	4505X	4.9	4.9	4.1	4.9	5.0	
WFF	4506X	4.0	4.0	4.0	4.0	4.9	
WFF	4507X	4.1	4.0	4.1	4.2	5.0	
WFF	4508X	3.0	4.0	4.0	3.9	4.9	small; PT
WFF	4509X	4.1	4.0	4.1	4.0	5.0	recheck hump; tight heart girth
WFF	4510X	4.9	4.1	4.1	4.2	4.0	
WFF	4511X	3.1	3.2	4.0	3.9	3.9	
WFF	4512X	4.9	4.1	4.9	4.9	5.0	
WFF	4513X	3.0	3.0	3.0	2.8	4.0	PT or CULL
WFF	4514X	3.1	3.9	3.9	3.9	5.0	
WFF	4515X						CULL
R2R	4613X	4.0	4.0	4.0	4.0	5.0	

CP	BULL	MASCU	THICK	MUSCLE	RATING	DISP	COMMENTS
R2R	4614X	3.0	2.9	3.0	2.9	5.0	tall, narrow
R2R	4615X	3.0	2.9	3.0	2.8	4.0	
R2R	4618X	4.0	4.0	4.1	4.1	4.9	
R2R	4619X	4.1	4.0	3.1	3.9	5.0	lame left rear; recheck
R2R	4628X	4.1	4.9	4.1	4.2	5.0	
R2R	4630X	3.0	3.1	3.1	3.0	5.0	
R2R	4638X	4.0	4.1	2.9	2.9	4.9	light rear; cow hock
R2R	4639X	3.0	3.0	2.9	2.9	4.0	dark red
R2R	4640X	3.9	3.9	3.9	3.2	4.9	
R2R	4650X	4.1	4.0	3.9	4.0	5.0	
R2R	4651X	3.0	3.1	3.1	3.1	4.9	
D-7	4726X	4.0	4.0	3.9	4.0	4.0	spot right eye
D-7	4728X	4.0	4.0	3.9	4.0	5.0	
D-7	4731X	3.1	3.9	3.1	3.0	5.1	dark red; cow hock
LLL	5901X	3.1	4.0	3.9	3.0	4.9	
LLL	5902X	4.1	4.9	4.1	4.2	5.0	cloudy left eye; right front twist; PT
LLL	5905X	4.1	4.1	4.0	3.9	5.0	recheck left rear
LLL	5907X	3.9	4.0	3.9	3.9	4.9	
LLL	5910X	4.9	5.0	5.0	5.0	5.1	
46							

CP	BULL	MASCU	THICK	MUSCLE	RATING	DISP	COMMENTS
D-J	2706X	3.9	3.9	3.9	3.8	4.9	spot right eye
D-J	2707X	3.1	3.0	2.9	3.0	3.9	
D-J	2708X	4.1	4.1	4.0	4.2	4.0	
D-J	2709X	5.0	3.0	3.9	3.1	3.9	tall
D-J	2710X	4.0	4.0	4.0	3.9	3.9	
D-J	2711X	2.9	2.9	2.9	2.2	3.9	
D-J	2712X	3.8	3.0	2.9	3.8	4.0	
D-J	2713X	5.0	4.1	3.8	4.9	5.0	
D-J	2714X	4.0	4.0	4.0	4.0	4.0	
D-J	2715X	4.8	4.1	3.0	4.2	4.9	
D-J	2716X	4.2	3.1	3.9	3.9	3.1	
D-J	2717X	4.0	4.2	3.9	4.1	5.0	DP
D-J	2718X	4.0	3.9	3.9	3.9	4.0	
D-J	2719X	2.9	2.9	3.0	2.9	4.1	
D-J	2720X	3.1	3.0	3.1	3.1	5.0	
D-J	2721X	4.0	3.9	3.1	3.9	3.1	
D-J	2722X	3.9	4.0	4.0	4.0	5.0	wide top
D-J	2723X	4.1	4.2	4.1	4.1	4.0	weak pasterns; rear toes spread; PT
D-J	2724X	3.1	4.0	4.0	4.0	4.1	
D-J	2725X	4.1	3.9	3.9	4.0	5.0	
D-J	2726X	3.9	3.1	3.0	3.0	4.0	spot left eye
D-J	2727X	4.1	4.1	4.0	4.1	4.1	
D-J	2729X	4.9	4.0	3.9	4.1	4.0	
D-J	2731X	3.1	2.9	2.9	2.9	4.9	spot & cloudy left eye
D-J	2732X	3.1	3.0	4.0	3.8	4.0	
D-J	2733X	3.0	3.0	3.9	4.0	2.1	
D-J	2734X	3.0	4.9	4.9	4.2	4.1	spot left eye

CP	BULL	MASCU	THICK	MUSCLE	RATING	DISP	COMMENTS
D-J	2735X	5.0	5.0	5.0	5.2	4.1	slope hip
D-J	2737X	3.1	3.0	3.0	3.0	3.0	
D-J	2738X	3.9	4.1	4.1	4.2	4.0	
D-J	2739X	3.1	4.0	4.0	3.2	4.9	spot left eye
D-J	2740X	3.0	3.2	4.0	3.9	4.9	
D-J	2743X	4.9	4.1	4.1	4.1	3.9	
D-J	2744X	4.0	3.1	3.0	3.0	4.9	
D-J	2745X	2.9	2.9	2.9	2.1	4.9	
D-J	2747X	3.9	3.8	3.2	3.1	3.9	
D-J	2750X	3.1	3.1	3.9	3.9	4.0	
D-J	2762X	4.1	4.9	4.1	5.0	5.0	
D-J	2764X	5.0	5.0	5.0	5.1	4.9	
39							